

Product Datasheet

Human Transient Receptor Potential Cation Channel Subfamily V, Member 1 (TRPV1) ELISA Kit (orb2810768)

Catalog Number	orb2810768
Category	Assays and Kits
Description	Human Transient Receptor Potential Cation Channel Subfamily V, Member 1 (TRPV1) ELISA Kit
Reactivity	Human
UniProt ID	Q8NER1
Application notes	Assay Type: Sandwich
Assay Type	Sandwich
Range	7.81-500pg/mL
Sensitivity	3.4pg/ml
Sample Types	Tissue homogenates;cell lysates and other biological fluids
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)