

Product Datasheet

Mouse Interferon Inducible T-Cell Alpha Chemoattractant (ITAC; CXCL11) ELISA Kit (orb2810791)

Catalog Number	orb2810791
Category	Assays and Kits
Description	Mouse Interferon Inducible T-Cell Alpha Chemoattractant (ITAC; CXCL11) ELISA Kit
Reactivity	Mouse
UniProt ID	Q9JHH5
Application notes	Assay Type: Sandwich
Assay Type	Sandwich
Range	31.2-2000pg/mL
Sensitivity	15.2pg/mL
Sample Types	Serum;plasma;tissue homogenates;cell lysates;cell culture supernates and other biological fluids
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)