

Product Datasheet

Human high sensitivity C-Reactive Protein (hs-CRP) ELISA Kit (orb2811115)

Catalog Number	orb2811115
Category	Assays and Kits
Description	Human high sensitivity C-Reactive Protein (hs-CRP) ELISA Kit
Reactivity	Human
UniProt ID	P02741
Application notes	Assay Type: Sandwich
Assay Type	Sandwich
Range	7.8-500pg/mL
Sensitivity	3.6pg/mL
Sample Types	Serum;plasma;tissue homogenates;cell lysates;urine;cerebrospinal fluid;cell culture supernates and other biological fluids
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)