

Product Datasheet

Human Sialate O-acetylsterases (SIAE) ELISA Kit (orb2811440)

Catalog Number	orb2811440
Category	Assays and Kits
Description	Human Sialate O-acetylsterases (SIAE) ELISA Kit
Reactivity	Human
UniProt ID	Q9HAT2
Application notes	Assay Type: Sandwich
Assay Type	Sandwich
Range	0.156-10ng/mL
Sensitivity	0.085ng/mL
Sample Types	Tissue homogenates;cell lysates and other biological fluids
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)