

Product Datasheet

Human Ubiquitin Specific Peptidase 14 (USP14) ELISA Kit (orb2811633)

Catalog Number	orb2811633
Category	Assays and Kits
Description	Human Ubiquitin Specific Peptidase 14 (USP14) ELISA Kit
Reactivity	Human
UniProt ID	P54578
Application notes	Assay Type: Sandwich
Assay Type	Sandwich
Range	0.312-20ng/mL
Sensitivity	0.15ng/mL
Sample Types	Tissue homogenates and other biological fluids
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)