

Product Datasheet

Human Cyclophilin of 40 kDa (CYP-40) ELISA Kit (orb2811972)

| | |
|--------------------------|---|
| Catalog Number | orb2811972 |
| Category | Assays and Kits |
| Description | Human Cyclophilin of 40 kDa (CYP-40) ELISA Kit |
| Reactivity | Human |
| UniProt ID | Q08752 |
| Application notes | Assay Type: Sandwich |
| Assay Type | Sandwich |
| Range | 0.312-20ng/mL |
| Sensitivity | 0.16ng/mL |
| Sample Types | Tissue homogenates;cell lysates;cell culture supernates and other biological fluids |
| Note | For research use only |
| Expiration Date | 6 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)