

Product Datasheet

Human Tissue Inhibitors Of Metalloproteinase 2 (TIMP-2) ELISA Kit (orb2812373)

Catalog Number	orb2812373
Category	Assays and Kits
Description	Human Tissue Inhibitors Of Metalloproteinase 2 (TIMP-2) ELISA Kit
Reactivity	Human
UniProt ID	P16035
Application notes	Assay Type: Sandwich
Assay Type	Sandwich
Range	0.31-20ng/mL
Sensitivity	0.15ng/mL
Sample Types	Serum;plasma;tissue homogenates and other biological fluids
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)