

Product Datasheet

Human Apelin 36 (AP36) ELISA Kit (orb2812804)

Catalog Number	orb2812804
Category	Assays and Kits
Description	Human Apelin 36 (AP36) ELISA Kit
Reactivity	Human
UniProt ID	Q9ULZ1
Application notes	Assay Type: Competitive
Assay Type	Competitive
Range	46.87-3000pg/mL
Sensitivity	12.58pg/mL
Sample Types	Serum;plasma and other biological fluids
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)