

Product Datasheet

Rat Trypsinogen Activation Peptide (TAP) ELISA Kit (orb2812932)

Catalog Number	orb2812932
Category	Assays and Kits
Description	Rat Trypsinogen Activation Peptide (TAP) ELISA Kit
Reactivity	Rat
UniProt ID	P00762
Application notes	Assay Type: Competitive
Assay Type	Competitive
Range	0.156-10ng/mL
Sensitivity	0.04ng/mL
Sample Types	Serum;plasma;tissue homogenates;cell lysates;cell culture supernates and other biological fluids
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)