

Product Datasheet

Rat Epithelial Neutrophil Activating Peptide 78 (ENA78) ELISA Kit (orb2812955)

Catalog Number	orb2812955
Category	Assays and Kits
Description	Rat Epithelial Neutrophil Activating Peptide 78 (ENA78) ELISA Kit
Reactivity	Rat
UniProt ID	P97885
Application notes	Assay Type: Sandwich
Assay Type	Sandwich
Range	15.62-1000pg/mL
Sensitivity	6.15pg/mL
Sample Types	Serum;plasma and other biological fluids
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)