

Product Datasheet

Human Ubiquitin Carboxyl Terminal Hydrolase L1 (Uch-L1) ELISA Kit (orb2813283)

Catalog Number	orb2813283
Category	Assays and Kits
Description	Human Ubiquitin Carboxyl Terminal Hydrolase L1 (Uch-L1) ELISA Kit
Reactivity	Human
UniProt ID	P09936
Application notes	Assay Type: Sandwich
Assay Type	Sandwich
Range	78-5000pg/mL
Sensitivity	37.06pg/mL
Sample Types	Tissue homogenates;cell lysates;cell culture supernates and other biological fluids
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)