

Product Datasheet

Human Adrenocorticotrophic Hormone (ACTH) ELISA Kit (orb2813443)

Catalog Number	orb2813443
Category	Assays and Kits
Description	Human Adrenocorticotrophic Hormone (ACTH) ELISA Kit
Reactivity	Human
UniProt ID	P01189
Application notes	Assay Type: Competitive
Assay Type	Competitive
Range	1.56-100pg/mL
Sensitivity	0.76pg/mL
Sample Types	Serum;plasma;tissue homogenates;cell lysates;cell culture supernates and other biological fluids
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)