

Product Datasheet

Mouse Cluster Of Differentiation 30 Ligand (CD30L) ELISA Kit (orb2813508)

Catalog Number	orb2813508
Category	Assays and Kits
Description	Mouse Cluster Of Differentiation 30 Ligand (CD30L) ELISA Kit
Reactivity	Mouse
UniProt ID	P32972
Application notes	Assay Type: Sandwich
Assay Type	Sandwich
Range	0.15-10ng/mL
Sensitivity	0.09ng/mL
Sample Types	Tissue homogenates;cell lysates and other biological fluids
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)