

Product Datasheet

Human CD5 Antigen Like Protein (CD5L) ELISA Kit (orb2813547)

Catalog Number	orb2813547
Category	Assays and Kits
Description	Human CD5 Antigen Like Protein (CD5L) ELISA Kit
Reactivity	Human
UniProt ID	O43866
Application notes	Assay Type: Sandwich
Assay Type	Sandwich
Range	0.78-50ng/mL
Sensitivity	0.29ng/mL
Sample Types	Serum;plasma;tissue homogenates;cell lysates and other biological fluids
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)