

Product Datasheet

Lymphotactin/XCL1 Protein, Mouse, Recombinant (orb2815918)

Catalog Number	orb2815918
Description	Lymphotactin/XCL1 Protein, Mouse, Recombinant is expressed in E. coli with Tag Free. The accession number is P47993.
Storage	-20°C
Tag	Tag Free
Note	For research use only
Application notes	Reconstitute the lyophilized protein in sterile deionized water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.
Protein Sequence	VGTEVLEESSCVNLQTQRLPVQKIKTYIIWEGAMRAVIFVTKRGLKICADPEAKWVKAAIKTV DGRASTRKNMAETVPTGAQRSTSTAITLTG
MW	10.1 kDa (Predicted)
Uniprot ID	P47993
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)