

Product Datasheet

CCL9 Protein, Mouse, Recombinant (Active) (orb2815975)

Catalog Number	orb2815975
Description	CCL9 Protein, Mouse, Recombinant (Active) is expressed in E. coli with Tag Free. The accession number is P51670.
Storage	-20°C
Tag	Tag Free
Note	For research use only
Application notes	Reconstitute the lyophilized protein in sterile deionized water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.
Protein Sequence	QITHATETKEVQSSLKAQQGLEIEMFHMGMFQDSSDCCLSYNSRIQCSRFIGYFPTSGGCTRP GIIFISKRGFQVCANPSDRRVQRCIERLEQNSQPRTYKQ
MW	11.6 kDa (Predicted)
Uniprot ID	P51670
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)