

Product Datasheet

CCL28 Protein, Human, Recombinant (Active) (orb2816361)

Catalog Number	orb2816361
Description	CCL28 Protein, Human, Recombinant (Active) is expressed in E. coli with Tag Free. The accession number is Q9NRJ3.
Storage	-20°C
Tag	Tag Free
Note	For research use only
Application notes	Reconstitute the lyophilized protein in sterile deionized water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.
Protein Sequence	SEAILPIASSCCTEVSHHISRLLERVNMCRIQRADGDCDLAAVILHVKRRRICVSPHNHTVK QWMKVQAAKNGKGNVCHRKKHHGKRNSNRAHQGKHETYGHKTPY
MW	12.4 kDa (Predicted)
Uniprot ID	Q9NRJ3
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)