

Product Datasheet

Human Holo TC (Holotranscobalamin) ELISA Kit (orb2816921)

Catalog Number	orb2816921
Description	Human Holo TC (Holotranscobalamin) ELISA Kit
Reactivity	Human
Range	0.625-40pmol/L(40pmol=1734ng)
Target	Holo TC
Storage	2 - 8 °C for 12 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, cell or tissue lysate, Other liquid samples
Assay Time	4 hours
Sensitivity	0.375pmol/L
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)