

## Product Datasheet

### n-Myc Antibody (orb2817055)

|                            |                                                                                                                                          |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>      | orb2817055                                                                                                                               |
| <b>Description</b>         | n-Myc Antibody                                                                                                                           |
| <b>Species/Host</b>        | Rabbit                                                                                                                                   |
| <b>Reactivity</b>          | Human, Mouse                                                                                                                             |
| <b>Conjugation</b>         | Unconjugated                                                                                                                             |
| <b>Tested Applications</b> | ICC, IP, WB                                                                                                                              |
| <b>Immunogen</b>           | between 100 and 150                                                                                                                      |
| <b>Target</b>              | n-Myc                                                                                                                                    |
| <b>Preservatives</b>       | Phosphate Buffered Saline (PBS) with 0.09% Sodium Azide, rAlbumin-Free                                                                   |
| <b>Concentration</b>       | 1000 µg/ml                                                                                                                               |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.  |
| <b>Note</b>                | For research use only                                                                                                                    |
| <b>Application notes</b>   | Format: Whole IgG                                                                                                                        |
| <b>Clonality</b>           | Recombinant                                                                                                                              |
| <b>Antibody Type</b>       | Primary Antibody                                                                                                                         |
| <b>Uniprot ID</b>          | <b>P04198</b>                                                                                                                            |
| <b>Dilution Range</b>      | WB: 1:1,000, IP: 6 µl/mg lysate, ICC: 1:100 to 1:500. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE cell sections |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                          |

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)