

Product Datasheet

phospho-Bcl-xL (Ser62) Rabbit pAb, PE-Cy5.5 conjugated (orb2822569)

Catalog Number	orb2822569
Category	Antibodies
Description	phospho-Bcl-xL (Ser62) Rabbit pAb, PE-Cy5.5 conjugated is a PE/Cy5.5 conjugated antibody targeting BCL2L1. This antibody is suitable for IF. It exhibits reactivity with Rat samples.
Target	BCL2L1
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	PE/Cy5.5
Reactivity	Rat
Predicted Reactivity	Bovine, Canine, Equine, Human, Mouse, Porcine, Sheep
Concentration	1mg/ml
Buffer/Preservatives	Preservative: 0.02% Proclin300, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4.
Purification	Affinity purified by Protein A
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Bcl-xL around the phosphorylation site of Ser62 AD(p-S)PA
Tested applications	IF

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Dilution range	IF=1:100-500
Antibody Type	Primary Antibody
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)