

Product Datasheet

Rat Glucoprotein 130 (gp130) ELISA Kit (orb2828178)

| | |
|------------------------|--|
| Catalog Number | orb2828178 |
| Category | Assays and Kits |
| Description | The Rat Glucoprotein 130 (gp130) ELISA Kit features a dynamic range of 3.12-200ng/mL and is suitable for Serum, plasma, tissue homogenates or other biological fluids. This product offers a sensitivity of 1.21ng/mL. |
| Reactivity | Rat |
| Range | 3.12-200ng/mL |
| Sensitivity | 1.21ng/mL |
| Sample Types | Serum, plasma, tissue homogenates or other biological fluids. |
| Storage | Store at 4°C & -20°C for 12 months. |
| Note | For research use only |
| Expiration Date | 6 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)