

Product Datasheet

Anti-human IL-1RAP / IL-1R3 (Nidanilimab Biosimilar) (orb2829108)

Catalog Number	orb2829108
Description	Anti-human IL-1RAP / IL-1R3 (Nidanilimab Biosimilar)
Reactivity	Human
Endotoxins	< 0.001 EU/ug
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, In vivo
Form/Appearance	Lyophilized
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Dissolve with sterile ddH2O
Isotype	Human IgG1
Clonality	Monoclonal
Antibody Type	Biosimilar Antibody
MW	144.48 kDa
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)