

## Product Datasheet

### human CD45RA Mouse mAb, PE-Cy5.5 conjugated (orb2830209)

<b>Catalog Number</b>	orb2830209
<b>Category</b>	Antibodies
<b>Description</b>	human CD45RA Mouse mAb, PE-Cy5.5 conjugated is a PE/Cy5.5 conjugated antibody targeting PTPRC. This antibody is suitable for FC. It exhibits reactivity with Human samples.
<b>Target</b>	PTPRC
<b>Clonality</b>	Monoclonal
<b>Species/Host</b>	Mouse
<b>Isotype</b>	Mouse IgG1, k
<b>Conjugation</b>	PE/Cy5.5
<b>Reactivity</b>	Human
<b>Buffer/Preservatives</b>	0.01M PBS (pH7.4), 0.5% rAlbumin, 0.02% Proclin300.
<b>Purification</b>	Affinity purified by Protein G
<b>Immunogen</b>	KLH conjugated synthetic peptide derived from human CD45RA
<b>Tested applications</b>	FC
<b>Dilution range</b>	Flow-Cyt=1ug/Test
<b>Antibody Type</b>	Primary Antibody
<b>Clone Number</b>	B67H3

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

<b>Storage</b>	Store at 2-8°C. Protect from light. Avoid freezing.
<b>Note</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.

---

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)