

Product Datasheet

human CD45RA Mouse mAb, PE-Cy5.5 conjugated (orb2830209)

Catalog Number	orb2830209
Category	Antibodies
Description	human CD45RA Mouse mAb, PE-Cy5.5 conjugated is a PE/Cy5.5 conjugated antibody targeting PTPRC. This antibody is suitable for FC. It exhibits reactivity with Human samples.
Target	PTPRC
Clonality	Monoclonal
Species/Host	Mouse
Isotype	Mouse IgG1, k
Conjugation	PE/Cy5.5
Reactivity	Human
Buffer/Preservatives	0.01M PBS (pH7.4), 0.5% rAlbumin, 0.02% Proclin300.
Purification	Affinity purified by Protein G
Immunogen	KLH conjugated synthetic peptide derived from human CD45RA
Tested applications	FC
Dilution range	Flow-Cyt=1ug/Test
Antibody Type	Primary Antibody
Clone Number	45F1

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)