

Product Datasheet

Recombinant Human CD118/LIFR Protein, N-His (orb2831910)

Catalog Number	orb2831910
Category	Proteins
Description	Recombinant Human CD118/LIFR Protein, N-His, expressed in E.coli, with a molecular weight of 18.19 kDa, and a purity of >90% as determined by SDS-PAGE. It was purified by AC, and has a N-His tag.
Target	Leukemia inhibitory factor receptor
Tag	N-His
Purity	>90% as determined by SDS-PAGE.
Purification	AC
Protein Sequence	45-184/1097
UniProt ID	P42702
MW	18.19 kDa
Application notes	Reconstitute in sterile water for a stock solution.
Expression System	E.coli
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)