

Product Datasheet

Recombinant Human CD147/BSG Protein, N-His (orb2831998)

| | |
|-----------------------------|--|
| Catalog Number | orb2831998 |
| Category | Proteins |
| Description | Recombinant Human CD147/BSG Protein, N-His, expressed in E.coli, with a molecular weight of 22.47 kDa, and a purity of >90% as determined by SDS-PAGE. It was purified by AC, and has a N-His tag. |
| Target | Basigin |
| Tag | N-His |
| Form/Appearance | Lyophilized |
| Buffer/Preservatives | Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol. |
| Reconstitution | Reconstitute in sterile water for a stock solution. |
| Purity | >90% as determined by SDS-PAGE. |
| Purification | AC |
| Protein Sequence | 140-323/385 |
| UniProt ID | P35613 |
| MW | 22.47 kDa |
| Expression System | E.coli |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

| | |
|------------------------|---|
| Storage | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. |
| Note | For research use only |
| Expiration Date | 6 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)