

Product Datasheet

Recombinant Mouse CD40/TNFRSF5 Protein, C-His (orb2834189)

Catalog Number	orb2834189
Category	Proteins
Description	Recombinant Mouse CD40/TNFRSF5 Protein, C-His, expressed in HEK293 cell, with a molecular weight of 22.39 kDa, and a purity of >90% as determined by SDS-PAGE. It was purified by AC, and has a C-His tag.
Target	Tumor necrosis factor receptor superfamily member 5
Tag	C-His
Form/Appearance	Lyophilized
Buffer/Preservatives	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Reconstitution	Reconstitute in sterile water for a stock solution.
Purity	>90% as determined by SDS-PAGE.
Purification	AC
Protein Sequence	1-193/289
UniProt ID	P27512
MW	22.39 kDa
Expression System	HEK293 cell
Biological Origin	Mouse

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage

Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note

For research use only

Expiration Date

6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)