

## Product Datasheet

### Recombinant Human GPR35 Protein, N-GST (orb2834553)

|                          |                                                                                                                                                                           |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>    | orb2834553                                                                                                                                                                |
| <b>Description</b>       | Recombinant Human GPR35 Protein, N-GST                                                                                                                                    |
| <b>Preservatives</b>     | Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.                                                                                |
| <b>Form/Appearance</b>   | Lyophilized                                                                                                                                                               |
| <b>Storage</b>           | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. |
| <b>Tag</b>               | N-GST                                                                                                                                                                     |
| <b>Note</b>              | For research use only                                                                                                                                                     |
| <b>Application notes</b> | Reconstitute in sterile water for a stock solution.                                                                                                                       |
| <b>Protein Sequence</b>  | 240-309/309                                                                                                                                                               |
| <b>Purity</b>            | > 90% as determined by SDS-PAGE.                                                                                                                                          |
| <b>MW</b>                | 34.46 kDa                                                                                                                                                                 |
| <b>Uniprot ID</b>        | <b>Q9HC97</b>                                                                                                                                                             |
| <b>Expiration Date</b>   | 6 months from date of receipt.                                                                                                                                            |

---

#### **Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### **Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)