

Product Datasheet

Recombinant Human IL36A/IL-1F6 Protein, N-His (orb2834739)

Catalog Number	orb2834739
Category	Proteins
Description	Recombinant Human IL36A/IL-1F6 Protein, N-His, expressed in E.coli, with a molecular weight of 19.85 kDa, and a purity of >90% as determined by SDS-PAGE. It was purified by AC, and has a N-His tag.
Target	Interleukin-36 alpha
Tag	N-His
Purity	>90% as determined by SDS-PAGE.
Purification	AC
Protein Sequence	1-158/158
UniProt ID	Q9UHA7
MW	19.85 kDa
Application notes	Reconstitute in sterile water for a stock solution.
Expression System	E.coli
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)