

Product Datasheet

Recombinant Human CD288/TLR8 Protein, N-His (orb2834802)

Catalog Number	orb2834802
Description	Recombinant Human CD288/TLR8 Protein, N-His
Preservatives	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Form/Appearance	Lyophilized
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Tag	N-His
Note	For research use only
Application notes	Reconstitute in sterile water for a stock solution.
Protein Sequence	27-140/1041
Purity	> 90% as determined by SDS-PAGE.
MW	15.26 kDa
Uniprot ID	Q9NR97
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)