

Product Datasheet

Recombinant Human CD273/PD-L2/PDCD1LG2 Protein, C-His (orb2834923)

Catalog Number	orb2834923
Category	Proteins
Description	Recombinant Human CD273/PD-L2/PDCD1LG2 Protein, C-His, expressed in HEK293 cell, with a molecular weight of 23.62 kDa, and a purity of >90% as determined by SDS-PAGE. It was purified by AC, and has a C-His tag.
Target	Programmed cell death 1 ligand 2
Tag	C-His
Form/Appearance	Lyophilized
Buffer/Preservatives	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Reconstitution	Reconstitute in sterile water for a stock solution.
Purity	>90% as determined by SDS-PAGE.
Purification	AC
Protein Sequence	1-220/273
UniProt ID	Q9BQ51
MW	23.62 kDa
Expression System	HEK293 cell
Biological Origin	Human

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Storage

Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note

For research use only

Expiration Date

6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)