

Product Datasheet

Recombinant Human IGFBP4 Protein, N-His (orb2836173)

Catalog Number	orb2836173
Category	Proteins
Description	Recombinant Human IGFBP4 Protein, N-His, expressed in E.coli, with a molecular weight of 28.28 kDa, and a purity of >90% as determined by SDS-PAGE. It was purified by AC, and has a N-His tag.
Target	Insulin-like growth factor-binding protein 4
Tag	N-His
Purity	>90% as determined by SDS-PAGE.
Purification	AC
Protein Sequence	22-258/258
UniProt ID	P22692
MW	28.28 kDa
Application notes	Reconstitute in sterile water for a stock solution.
Expression System	E.coli
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)