

Product Datasheet

Recombinant human CD30 protein, (HEK293) C-His&Avi (orb2837206)

Catalog Number	orb2837206
Description	Recombinant human CD30 protein, (HEK293) C-His&Avi
Endotoxins	Not analyzed
Target	Tumor necrosis factor receptor superfamily member 8
Preservatives	PBS (pH=7.4)
Form/Appearance	Liquid
Concentration	> 1mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-10×His & C-Avi
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	19-379/595
Purity	> 90% as determined by SDS-PAGE
MW	42.1 kDa
Uniprot ID	P28908
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)