

Product Datasheet

phospho-IKK alpha (Tyr463) Rabbit pAb, Pacific Blue conjugated (orb2841075)

Catalog Number	orb2841075
Description	phospho-IKK alpha (Tyr463) Rabbit pAb, Pacific Blue conjugated
Species/Host	Rabbit
Reactivity	Human
Conjugation	Pacific Blue
Tested Applications	FC, ICC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human IKK alpha around the phosphorylation site of Tyr463 NL(p-T)KM
Target	CHUK
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Antibody Type	Primary Antibody
Dilution Range	ICC/IF=1:50-200, Flow-Cyt=1ug/test
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+441223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)