

## Product Datasheet

### PADI4 Rabbit pAb, PerCP-Cy7 conjugated (orb2850507)

|                             |                                                                                                                                                                         |
|-----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>       | orb2850507                                                                                                                                                              |
| <b>Category</b>             | Antibodies                                                                                                                                                              |
| <b>Description</b>          | PADI4 Rabbit pAb, PerCP-Cy7 conjugated is a PerCP/Cy7 conjugated antibody targeting PADI4. This antibody is suitable for IF. It exhibits reactivity with Human samples. |
| <b>Target</b>               | PADI4                                                                                                                                                                   |
| <b>Clonality</b>            | Polyclonal                                                                                                                                                              |
| <b>Species/Host</b>         | Rabbit                                                                                                                                                                  |
| <b>Isotype</b>              | IgG                                                                                                                                                                     |
| <b>Conjugation</b>          | PerCP/Cy7                                                                                                                                                               |
| <b>Reactivity</b>           | Human                                                                                                                                                                   |
| <b>Concentration</b>        | 1mg/ml                                                                                                                                                                  |
| <b>Buffer/Preservatives</b> | 0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.                                                                                                  |
| <b>Purification</b>         | Affinity purified by Protein A                                                                                                                                          |
| <b>Immunogen</b>            | KLH conjugated synthetic peptide derived from human PADI4 (401-500/663aa)                                                                                               |
| <b>Tested applications</b>  | IF                                                                                                                                                                      |
| <b>Dilution range</b>       | IF=1:100-500                                                                                                                                                            |
| <b>Antibody Type</b>        | Primary Antibody                                                                                                                                                        |

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Storage**

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

**Note**

For research use only

---

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)