

Product Datasheet

phospho-MERTK (Tyr749 + Tyr753 + Tyr754) Rabbit pAb, IRDye 800CW conjugated (orb2855299)

Catalog Number	orb2855299
Category	Antibodies
Description	phospho-MERTK (Tyr749 + Tyr753 + Tyr754) Rabbit pAb, IRDye 800CW conjugated is a IRDye800 conjugated antibody targeting MERTK. This antibody is suitable for WB. It exhibits reactivity with Human samples.
Target	MERTK
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	IRDye800
Reactivity	Human
Predicted Reactivity	Mouse, Rat
Concentration	1mg/ml
Buffer/Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Purification	Affinity purified by Protein A
Immunogen	KLH conjugated synthesised phosphopeptide derived from human MERTK around the phosphorylation site of Tyr749 + Tyr753 + Tyr754 KI(p-Y)SGD(p-Y)(p-Y)RQ

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Tested applications	WB
Dilution range	WB=1000-10000
Antibody Type	Primary Antibody
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)