

Product Datasheet

phospho-PDCD4 (Ser457) Rabbit pAb, Pacific Blue conjugated (orb2862655)

Catalog Number	orb2862655
Category	Antibodies
Description	phospho-PDCD4 (Ser457) Rabbit pAb, Pacific Blue conjugated is a Pacific Blue conjugated antibody targeting PDCD4. This antibody is suitable for ICC. It exhibits reactivity with Human, Mouse samples.
Target	PDCD4
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	Pacific Blue
Reactivity	Human, Mouse
Predicted Reactivity	Equine, Gallus, Porcine, Rabbit, Rat, Sheep, Zebrafish
Concentration	1mg/ml
Buffer/Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Purification	Affinity purified by Protein A
Immunogen	KLH conjugated synthesised phosphopeptide derived from human PDCD4 around the phosphorylation site of Ser457 FV(p-S)EG
Tested applications	ICC

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Dilution range	ICC/IF=1:50-200
Antibody Type	Primary Antibody
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)