

Product Datasheet

MALT1 Rabbit pAb, AE conjugated (orb2869612)

Catalog Number	orb2869612
Category	Antibodies
Description	MALT1 Rabbit pAb, AE conjugated
Target	MALT1
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	AE
Predicted Reactivity	Equine, Gallus, Human, Mouse, Rat, Sheep
Concentration	1mg/ml
Buffer/Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Purification	Affinity purified by Protein A
Immunogen	KLH conjugated synthetic peptide derived from human MALT1 (631-730/824aa)
Tested applications	ELISA
Dilution range	ELISA=1:500-5000
Antibody Type	Primary Antibody
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558