

Product Datasheet

phospho-MYPT1 (Thr696) Rabbit pAb, BF700 conjugated (orb2884983)

Catalog Number	orb2884983
Category	Antibodies
Description	phospho-MYPT1 (Thr696) Rabbit pAb, BF700 conjugated is a BF700 conjugated antibody targeting PPP1R12A. This antibody is suitable for FC, ICC. It exhibits reactivity with Human samples.
Target	PPP1R12A
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	BF700
Reactivity	Human
Predicted Reactivity	Bovine, Canine, Equine, Gallus, Mouse, Rabbit, Rat
Concentration	1mg/ml
Buffer/Preservatives	Preservative: 0.02% Proclin300, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4.
Purification	Affinity purified by Protein A
Immunogen	KLH conjugated synthesised phosphopeptide derived from human MYPT1 around the phosphorylation site of Thr696 RS(p-T)QG
Tested applications	FC, ICC

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Dilution range	ICC/IF=1:50-200, Flow-Cyt=1ug/test
Antibody Type	Primary Antibody
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)