

## Product Datasheet

# HA Tag Recombinant Rabbit mAb, Biotin conjugated (orb2891835)

<b>Catalog Number</b>	orb2891835
<b>Category</b>	Antibodies
<b>Description</b>	The HA Tag Recombinant Rabbit mAb, Biotin conjugated is a Biotin conjugated Recombinant antibody suitable for WB. This product reacts with multiple species.
<b>Clonality</b>	Recombinant
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	Biotin
<b>Reactivity</b>	All
<b>Purification</b>	Affinity purified by Protein A
<b>Immunogen</b>	KLH conjugated HA Tag
<b>Tested applications</b>	FC, ICC, IF, IHC-Fr, IHC-P, WB
<b>Dilution range</b>	WB=1:5000-50000, IHC-P=1:200-1000, IHC-F=1:200-1000, IF=1:200-1000, Flow-Cyt=1ug/Test, ICC/IF=1:100-500
<b>Antibody Type</b>	Primary Antibody
<b>Storage</b>	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)