

Product Datasheet

GST tag Recombinant Rabbit mAb, HRP conjugated (orb2891858)

Catalog Number	orb2891858
Category	Antibodies
Description	The GST tag Recombinant Rabbit mAb, HRP conjugated is an HRP conjugated Recombinant antibody suitable for WB. This product reacts with multiple species.
Clonality	Recombinant
Species/Host	Rabbit
Conjugation	HRP
Reactivity	All
Buffer/Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Purification	Affinity purified by Protein A
Immunogen	Recombinant GST Protein
Tested applications	IHC-Fr, IHC-P, WB
Dilution range	WB=1:500-2000, IHC-P=1:200-1000, IHC-F=1:200-1000
Antibody Type	Primary Antibody
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)