

Product Datasheet

Mouse Calcitonin Gene Related Peptide (CGRP) ELISA Kit (orb2892998)

Catalog Number	orb2892998
Category	Assays and Kits
Description	The Mouse Calcitonin Gene Related Peptide (CGRP) ELISA Kit features a dynamic range of 15.6-1000pg/mL and is suitable for Serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. This product offers a sensitivity of 4.22pg/mL.
Reactivity	Mouse
Range	15.6-1000pg/mL
Sensitivity	4.22pg/mL
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Storage	Store at 4°C & -20°C for 12 months.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)