

Product Datasheet

RPA32, Phospho (S33) Antibody (orb2893114)

Catalog Number	orb2893114
Description	RPA32, Phospho (S33) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC
Immunogen	surrounding Serine 33
Target	RPA32, Phospho (S33)
Preservatives	Phosphate Buffered Saline (PBS) with 0.1% rAlbumin and 0.09% Sodium Azide
Concentration	100 µg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Recombinant
Antibody Type	Primary Antibody
Uniprot ID	P15927
Dilution Range	ICC: 1:100 to 1:500. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE cell sections
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)