

## Product Datasheet

### mouse IL-1b(10G11) Antibody (orb2899201)

<b>Catalog Number</b>	orb2899201
<b>Category</b>	Antibodies
<b>Description</b>	The 10G11 monoclonal antibody is a mouse anti-mouse reagent designed for the specific detection and neutralization of interleukin-1 beta (IL-1 $\beta$ ). It is widely used in both in vitro assays and in vivo models to study inflammation, autoimmune diseases, and the IL-1 signaling pathway.
<b>Target</b>	IL Receptor
<b>Conjugation</b>	Unconjugated
<b>Antibody Type</b>	Inhibitor Antibody
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)