

Product Datasheet

mouse IL-1b Antibody (orb2899201)

Catalog Number	orb2899201
Category	Antibodies
Description	The B32I33 monoclonal antibody is a mouse anti-mouse reagent designed for the specific detection and neutralization of interleukin-1 beta (IL-1 β). It is widely used in both in vitro assays and in vivo models to study inflammation, autoimmune diseases, and the IL-1 signaling pathway.
Target	IL Receptor
Conjugation	Unconjugated
Antibody Type	Inhibitor Antibody
Clone Number	B32I33
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)