

Product Datasheet

Recombinant Human Furin(FURIN) (orb2922873)

Catalog Number	orb2922873
Description	Recombinant Human Furin(FURIN)
Preservatives	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Form/Appearance	Lyophilized powder
Storage	Storage Condition: Store at -20°C upon receipt, aliquoting is necessary for mutiple use. Avoid repeated freeze-thaw cycles. Shelf Life: The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. Notes: Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week
Note	For research use only
Protein Sequence	DVYQEPTDPKFPQQWYLSGVTQRDLNVKAAWAQGYTGHGIVVSILDDGIEKNHPDLAGNY DPGASFVDNDQDPDPQPRYTQMNDNRHGTRCAGEVAAVANNGVCGVGVAYNARIGGV ML DGEVTDAVEARSLGLNPNHIHIYSASWGPEDDGKTVDPARLAEAAFFRGVSQGRGGLGS IFVWASGNGGREHDSNCNCDGYTNSIYTLSSSATQFGNVPWYSEACSSTLATTYSSGNQN EKQIVTTDLRQKCTESHTGTSASAPLAAGIALTLANKNLTWRDMQHLVVQTSKPAHLN ANDWATNGVGRKVSHSYGYLLDAGAMVALAQNWTTVAPQRKCIIDILTEPKDIGKRLEV RKTVTACLGEPNHITRLEHAQARLTLSYNRRGDLAIHLVSPMGTRSTLLAARPHDYSADG FNDWAFMTTHSWDEDPSGEWVLEIENTSEANNYGTLTKFTLVLYGTAEGLPVPPESSGC KTLTSSQACVVCEEGFSLHQKSCVQHCPPGFAPQVLDTHYSTENDVETIRASVCAPCHAS CATCQGPALTDCLSCPSHASLDPVEQTCSRQSQSSRESPPQQPPRLPPEVEAGQRLRAG LLPSHLPEVVAGLSCAFIVLVFVTVFLVLQLRSGFSFRGVKYYTMDRGLISYKGLPPEAW QEECPSDSEEDGREGERTAFIKDQSAL
Uniprot ID	P09958
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)