

Product Datasheet

Bovine Genomic DNA (orb2937600)

Catalog Number	orb2937600
Category	Small Molecules
Description	Bovine genomic DNA is recommended for genomic analysis e.g. as template in positive control PCRs like for the MeatDetect qPCR Kit Pork+Horse+Beef (Halal), #PCR-708. The DNA was purified from bovine muscle meat. The DNA fragments have a size of up to 50 kb.
Form/Appearance	liquid
Concentration	200 ng/μl
Storage	store at -20 °C, avoid freeze/thaw cycles, aliquoting of DNA samples is recommended
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)