

Product Datasheet

Human Interleukin 12,IL-12/P40 ELISA KIT (orb2937720)

Catalog Number	orb2937720
Category	Assays and Kits
Description	This is the Human Interleukin 12 kit. It can be used with serum, plasma, tissue homogenates samples. The detection range spans 7.8 pg/mL to 500 pg/mL, offering a sensitivity of 3.9 pg/mL.
Reactivity	Human
UniProt ID	P29460
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	7.8pg/mL-500pg/mL
Sensitivity	3.9pg/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)