

Product Datasheet

Rabbit Anti-Sheep IgG (H+L) Secondary Antibody (FITC) (orb2938645)

Catalog Number	orb2938645
Category	Antibodies
Description	Rabbit Anti-Sheep IgG (H+L) Secondary Antibody (FITC)
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	FITC
Reactivity	Sheep
Form/Appearance	liquid
Buffer/Preservatives	5mg/ml rAlbumin (IgG-Free, Protease-Free)
Tested applications	FC, IF
Dilution range	1:20-1:100 for most applications
Antibody Type	Secondary Antibody
Storage	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol,pH7.6. -20°C for 24 months.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558