

Product Datasheet

TLR7 agonist 30 (orb2944677)

Catalog Number	orb2944677
Category	Small Molecules
Description	TLR7 agonist 30 (compound 2) is a cleavable linker-drug conjugate for ADC development, derived from TLR7 agonist 1. It demonstrates robust anti-tumor efficacy in research models, applicable for studying immune-stimulating antibody conjugates in both in vitro and in vivo oncology research.
Target	TLR, Drug-Linker Conjugates for ADC
MW	1529.77
Formula	$C_{72}H_{116}N_{14}O_{22}$
Storage	-20°C
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)